

Dornier **HOOVER**

FANS Ureteral Access Sheath




flexible sheath tip



control the negative pressure by adjusting the opening of vent



urological sampler set



hydrophilic coating on sheath and dilator

order information

Product Code	Product Description	Qty
IVX-NS-1040	Hoover Negative Pressure Ureteral Access Sheath, 10Fr x 40 cm	1/box
IVX-NS-1050	Hoover Negative Pressure Ureteral Access Sheath, 10Fr x 50 cm	1/box
IVX-NS-1140	Hoover Negative Pressure Ureteral Access Sheath, 11Fr x 40 cm	1/box
IVX-NS-1150	Hoover Negative Pressure Ureteral Access Sheath, 11Fr x 50 cm	1/box
IVX-NS-1240	Hoover Negative Pressure Ureteral Access Sheath, 12Fr x 40 cm	1/box
IVX-NS-1250	Hoover Negative Pressure Ureteral Access Sheath, 12Fr x 50 cm	1/box
IVX-NB-150	Hoover Urological Sampler Set	1 unit



Dornier MedTech