

Dornier
AXIS LEO

Clean & Clear, Every Time

Single-use digital flexible ureterorenoscope



9.6

Take
your
pick

7.5

AVAILABLE IN 2 SIZES



Dornier MedTech



If size matters,

make your choice without compromise on performance.

Available in 2 sizes,
7.5 Fr & 9.6 Fr



Why AXIS LEO

Brilliant imaging capabilities

- No reprocessing
- No repairs
- Readily available
- Always clean
- Cost-effective





| Video Processor

- Multiple connectivity options (DVI, S-Video, CVBS) allow compatibility with most monitor types for easy integration into any OR
- User-friendly menu navigation
- Integrated USB data storage option

1 270°
bidirectional
deflection

2 Passive angulation
supports access to
difficult anatomical
regions

3 Torque-stable
shaft assists
controlled
navigation

4 Luer-lock connector
(included) allows
the irrigation's
connection to the
working channel
access

5 Controls for white
balance and
image recording
positioned right at
your fingertips

3

2

AXIS LEO Flexible Video Ureterorenoscope

	9.6 Fr	7.5 Fr
Order information	K1046369 9.6 Fr URS EU deflection	K1048308 7.5 Fr URS EU deflection
	K1046370 9.6 Fr URS US deflection	K1048309 7.5 Fr URS US deflection
Distal tip	9.6 Fr / 3.2 mm	7.5 Fr / 2.5 mm
Working channel	3.6 Fr / 1.2 mm	
Working length	650 mm	
Deflection	270° up / down	
Illumination	Optical Fiber	Dual LED
Field of view	120°	
Depth of field	3 - 50 mm	
CMOS pixel	400 * 400	

AXIS LEO Video Processor

Order information	K1046368
Dimensions	(W) 24 cm x (L) 30 cm x (H) 6 cm
Weight	0.75 kg
Connections	USB 2.0, DVI, S-Video, CVBS



Scan the QR code for our
global office locations

www.dornier.com/locations